

BROADBAND PLASMA LIGHT SOURCE XWS-30








Operation and safety manual

Contents

| | |
|---------------------------------------------------------|----|
| 1. Content of the system..... | 3 |
| 2. Safety precautions | 4 |
| 3. Preparation for operation | 5 |
| 4. LED indication and system control panel..... | 6 |
| Buttons:..... | 6 |
| Indicators: | 6 |
| 5. System start and normal operational check-list | 7 |
| 6. System remote control via XWS-Monitor | 9 |
| 7. Control via RS terminal | 13 |
| 8. Specifications and facility requirements..... | 15 |
| 9. XWS-30 source dimensions..... | 16 |
| 10. Software download | 16 |
| 11. ISTEQ contacts | 17 |
| 12. Photonics4Work contacts | 17 |

1. Content of the system

ISTEQ broadband plasma light source XWS-30 system includes:

| | | | |
|----------------------------------|--|--------------------------------------------------------------------------------------|--|
| XWS-30 Optical head, x1 | |  | |
| Power supply unit (PSU) , x1 | |  | |
| PSU power cable, x1 | |  | |
| Interlock, x1 | |  | |
| RS-485/USB converter + cable, x1 | |  | |
| PC - Converter USB cable, x1 | |  | |
| Purge fittings, x2 | |  | |

After unpacking the system please check that there is no damage. If there is - contact the supplier.

2. Safety precautions

The XWS-30 laser plasma source uses a high-power laser light which is delivered by fiber and focused on plasma inside the light source head unit.

The plasma light source is very bright and can potentially cause damage to the eyes and skin if there is direct contact.



Exploitations restrictions:

- **DO NOT** open the light source head
- **DO NOT** touch the output window and the lamp bulb inside the light source head unit
- Operate the source in faculties with ambient temperature below 30°C
- Do not restrict air convection of the Optical head. Minimum distant to the nearest obstacles:
 - From source head back: 20cm
 - From source front panel: 10cm
- Light source head module gets HOT (up to 60°C under normal conditions during the first 45 minutes of continuous operation)

For your own protection and safety:

- Wear UV-protection glasses
- Take necessary precautions to avoid UV exposure
- Limit exposure to UV-generated ozone

3. Preparation for operation

1. Carefully unpack the delivered system and prepare for the operation:



2. Remove the plastic/metal cap from the optical head output window/FCU:



3. Connect the interlock (1), RS-485 cable (2) and the power cable (3) to the XWS-30 system:



4. Connect the power cable to 220V plug

5. Connect the RS-485/USB converter to your PC/Laptop via USB cable.




6. The system is ready for operation

4. LED indication and system control panel

LED indicators and control buttons are located on the side of the Optical head.

Buttons:

- ON/OFF button is used to switch ON and OFF the source
- START/STOP button is used to START and STOP the plasma



Indicators:

The indication system consists of five LEDs, two of those are integrated into the controlled buttons.



Turn ON the system pushing ON/OFF button one time, check the indication:


| READY LED | |
|-------------|---------------------------------------------------------------------------|
| No light | The system is not ready for the operation |
| GREEN light | The system is ready for operation |
| TEMP LED | |
| No light | Temperature control system is switched OFF and there is no internal error |

| | |
|------------------------------|----------------------------------------------------------------------------------------------------------------|
| RED light | Temperature control system is switched OFF and there is an internal error (In that case contact your supplier) |
| GREEN blinking light | Temperature control system is active and it is adjusting the laser temperature |
| GREEN light | Temperature of the laser is stabilized |
| LASER LED | |
| No light | The laser is switched OFF |
| RED light | The laser is switched OFF or there is an internal error |
| GREEN light | The laser is switched ON |
| ON/OFF LED/button | |
| No light | The system is powered OFF |
| GREEN light | The system is power ON |
| START/STOP LED/button | |
| No light | The plasma is switched OFF |
| GREEN blinking | Plasma ignition is ON |
| GREEN light | Plasma is ON |

5. System start and normal operational check-list

When the system is prepared for the operation do the following:

| | |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|
| 1. Connect the interlock (1), RS-485 cable(2) [Optional] and the power cable (3) to the XWS-30 system: |  |
| 2. Connect the power cable to 220V plug | |
| 3. Push ON/OFF button one time. Check that: <ul style="list-style-type: none"> - ON/OFF LED is green - READY LED is green - TEMP LED is green |  |


| | |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------|
| <p>4. Check there is no protection cap on the output flange:</p> |  |
| <p>5. Start the plasma pushing START/STOP button. Check that:</p> <ul style="list-style-type: none">- Plasma light is produced- TEMP LED is Green- There is no RED LED on the indication panel | |

6. System remote control via XWS-Monitor


XWS-Monitor is a special software which allows to:

- Remotely control the system
- Turn the plasma ON and OFF
- Check the system parameters online
- Log internal system parameters in a log file


To get the latest version of the «XWS-Monitor» software please contract our supplier or ISTEQ directly.

| | |
|-----------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| <p>1. Connect the interlock (1), RS-485 cable (2) and the power cable (3) to the XWS-30 system:</p> |  |
| <p>2. Connect the power cable to 220V plug.</p> | |
| <p>3. Contact your supplier to get the latest version of “XWS-30 Monitor” software.</p> | |
| <p>4. Turn ON the XWS-30 source and start “XWS-30 Monitor” software.</p> | |


- Click RESCAN button, choose the available COM port. Click CONNECT button. If the system is now connected you will send the system Serial number next to S/N:



- Go to SETTINGS tab. Choose the measurement time, which can be anything down to 200ms. Click APPLY




- TABLE tab. Here you can monitor the system parameters in real time:




8. Plots tab

Here you can see in real time the following graphs:


- Laser temperature, °C
- Laser current, A
- Plasma Photodiode level, mV
- Head case temperature, °C
- TEC (laser cooling system temperature), °C



9. Use START, STOP and RESTART buttons to start, stop the plasma or to restart the system:



10. Click LOG if you want the system parameters to be automatically saved in .txt file in the root folder of XWS-Monitor




List of system parameters can be seen via XWS-30 Monitor:

| Parameter | Meaning | Normal value |
|--------------------------|----------------------------------------------|---------------------------------------------------------------|
| Serial # | Serial number of the system | XXXX |
| Firmware version | Firmware version of the PSU control PCB | 1.11 or higher |
| Uptime | Total Uptime of the system, hh:mm | NA |
| Power On Time | Time of the current system session, hh:mm:ss | NA |
| Laser temperature | Temperature of the laser, °C | 25±0.1°C |
| Head temperature | Temperature of the Optical Head | 15°C < Th < 80°C |
| TEC current | Laser Cooling system: TEC element current | <8A |
| TEC voltage | Laser Cooling system: TEC element voltage | <15V |
| FAN Speed | Laser Cooing system: FAN rotation speed | <100% |
| Laser current | Current of the drive laser | 10±0.1A or 11±0.1A (depending on the system configuration) |
| Laser voltage | Voltage of the drive laser | <5V |

7. Control via RS terminal

That is possible to control the XWS-30 light source using RS-485 protocol.

Next steps to be taken:

| | |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| <p>1. Connect your PC/Laptop to XWS-30 system using RS-485/USB connector:</p> | |
| <p>2. It is up to the customer – to decide which software to be used for the RS communication As an option – Free Software <u>“Terminal”</u></p> |  |
| <p>3. Serial port settings: 115200 8-N-1</p> | |
| <p>4. Choose the COM-port and click <i>Connect</i></p> | |
| <p>5. Send commands from the list below to receive the information about the system.</p> <p>NOTE: That is also possible to START and PLASMA remotely using RS protocol. See the commands description below</p> | |

Command list and comments are shown in the table below:


| Command | Response |
|--------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| STATUS | STATUS=x where 'x' is system mode: x = 0 – IDLE (ready for operations) x = 1 – STARTING (Laser is starting) x = 2 – IGNITION (Plasma is triggered) x = 3 – PLASMA ON (Plasma is turned ON) x = 4 – ERROR |
| REPORT | The system gives the list of source internal parameters |
| ERROR | ERROR=xxxxx Allows to know error flags if STATUS == 4 0000000 means NO ERROR |
| SERIAL | SERIAL=xxxx where 'xxxx' is the system serial number shown as 4-digit number |
| FIRMWARE | FIRMWARE=x.xx where 'x.xx' is controller firmware version |
| UPTIME | UPTIME=xxxx.x where 'xxxx.x' is total laser/plasma active time in hours |
| PWRONTIME | PWRONTIME=xxxx days xx hours xx minutes xx seconds where 'xxxx days xx hours xx minutes xx seconds' is time since the system was switched on |
| LASER_TEMP | LASER_TEMP=xx.xxxx where 'xx.xxxx' is laser module temperature , °C |
| HEAD_TEMP | HEAD_TEMP=xx.xxxx where 'xx.xxxx' is optical head temperature, °C |
| TEC_CUR1 | TEC_CUR1=xx.xxxx where 'xx.xxxx' is TEC channel current in Amps |
| LASER_CUR | LASER_CUR=xx.xxxx where 'xx.xxxx' is laser current is Amps |
| LASER_VOL | LASER_VOL=xx.xxxx where 'xx.xxxx' is laser voltage in Volts |
| SUPPLY_VOL | SUPPLY_VOL=xx.xxxx where 'xx.xxxx' is system supply voltage in Volts |
| TURN_ON | TURN_ON=[WAIT OK ERROR] Turns plasma on & returns result of operation |
| TURN_OFF | TURN_OFF=[OK ERROR] Turns plasma off & returns result of operation |
| LASER_STAT | LASER=[ON OFF ERROR] |
| PLASMA_STAT | PLASMA=[ON OFF ERROR] |
| RESTART | Restarts the system |
| POWEROFF | Powers OFF the system |

8. Specifications and facility requirements

| XWS-30 performance | |
|----------------------------------|---------------------------------------------------------------------------------|
| Spectral range | 190 to 2500 nm (UV configuration), 250 to 2500 nm (Ozone-free configuration) |
| Spectral brightness (450-500 nm) | Up to 40mW/(mm ² ×sr×nm) |
| Emitting body source dimensions | 100×250μm |
| Lifetime | 10,000 hours |
| Temporal and spatial stability | STD < 0.15% |
| Main parameters | |
| Light source dimensions | 120x150x144.5mm |
| Output NA (by default) | 0.4 |
| External optic interface | C-mount |
| Pumping laser wavelength | NIR |
| Optional configurations: | |
| Source spectrum | UV or Ozone free |
| Light output | Free Space or Fiber coupled |
| Additional | |
| External control | COM-port (RS-485) |
| Interlock | LEMO FGG |
| Electricity requirements | |
| Voltage | 110-220V±10% |
| Frequency | 50Hz |

9. XWS-30 source dimensions

Linear dimensions: 120x150x144.5mm



10. Software download

The XWS-30 software is available as download from

www.photonics4work.eu/Download/XwS30/XWS-30.exe

Latest version of this manual

www.photonics4work.eu/Download/XwS30/manual.pdf

11. ISTEQ contacts

If you have any questions regarding the XWS source – feel free to contact your supplier or ISTEQ directly:

Address: The Netherlands, Eindhoven, High Tech Campus 9, 5656AE

Website: www.isteq.nl

Phones: +3140-230-42-42 (Office), +316-2525-7382 (Mobile)

E-mail: info@isteq.nl

12. Photonics4Work contacts

For any commercial information or if you have any questions regarding the XWS source contact:

Address: The Netherlands, Zevenaar, 6903 PZ, Mercurion 28 A

Website: www.photonics4work.eu

Phones: +31 31 – 316 340804 (Office), +31 6 – 22 40 60 27 (mobile)

E-mail: contact@photonics4work.eu